

Data Processing in MNova NMR


Below is a series of video tutorials covering different topics of MNova NMR processing. It is advisable to watch them in a sequence as they refer to previously discussed tools.

Episode 1. Which MNova do I need?

Removed.

Episode 2. Creating MNova Document

This video covers basic toolbars and menus of the MNova window. We create a new document with a title page and save it.




Summary

- Make a title page with the basic sample info for each document
 - Save with the original NMR data
 - Save the MNova document frequently as you progress through your project
-

Episode 3. Load data and move around with zoom tools

This episode demonstrates how to load NMR data into MNova document and explore the NMR spectrum using zoom tools.




Summary

- Load NMR data through drag-and-drop or Data Browser
- Update the spectrum Comment field and adjust its fonts size for display
- Add a title to the page itself to be able to find it easier in the Pages view
- Save frequently! (Cmd-S or Ctrl-S)
- Zoom intensity with a mouse wheel
- Zoom with 'Z', unzoom with Shift-'Z'; repeat a keystroke to switch the zooming direction
- Scale up with 'H' and 'F'. Pan with P.
- Navigate a spectrum with Shift-Z-Z-Z mode: click to zoom out, drag a box to zoom in
- Hit ESC on a keyboard to exit a zoom mode (or any other MNova tool)


Episode 4. Manual Phase and Baseline Correction

This episode demonstrates a brief workflow for manual phase correction in a full version of MNova NMR.



1D NMR Processing: Reference Deconvolution

This episode introduces a way to improve peak shapes in the spectrum if they were distorted by imperfect shimming. The tutorial focuses on basic concepts and a practical workflow. For more background on reference deconvolution, see this [blog post](#)



Work in progress! Future episodes are coming...

From:

<https://pydio.campus.nd.edu/docs/> - **NMR**

Permanent link:

https://pydio.campus.nd.edu/docs/doku.php?id=mnova_data_processing

Last update: **2025/05/09 19:21**

